




Olink®
Accelerating proteomics together

Olink Explore HT

Capture true biological
insights with proven
specificity. At any scale.



Olink®
Accelerating proteomics together



Olink Explore HT

5,400+ protein biomarkers

Increased library

Simplified workflow

Higher throughput

Automated data processing



Olink®
Accelerating proteomics together

Unprecedented power

Explore thousands of proteins with high specificity, minimal sample consumption, and user-friendly data tools for projects of any size



Proteins

5,400+

with proven specificity



Sample

2 μ L

Plasma, serum & more



Throughput

4x

Increase



Dynamic range

fg-mg/mL

10 logs



Runs on

NGS

Automated workflow



Olink®
Accelerating proteomics together

Measure proteins with confidence, obtaining reproducible results

Mean intra-CV
(within run)

<11%

Mean inter-CV
(between runs)

<9%

Per-assay correlation coefficient (R)
between Olink Explore 3072 and Olink Explore HT

0.88 median



Olink®
Accelerating proteomics together



Components

6x fewer

saves time and resources

Boxes

10x fewer

benefits the environment and saves
freezer space



Olink®
Accelerating proteomics together

Significant workflow cost savings

~20%

reduction in workflow
cost per sample

4x

reduction in
labor cost

4x

reduction in
consumables cost

Compared to Olink Explore 3072. The exact cost reduction varies depending on the labor cost.



Olink®
Accelerating proteomics together

80% more data at same sequencing cost/sample



Runs on Illumina® NovaSeq 6000



344 biological samples per run (two S4 flow cells)

Compared to Olink Explore 3072



Olink®
Accelerating proteomics together

Re-imagine your workflow, not your lab

No new instrumentation required




dragonfly® discovery
SPT Labtech



F.A.S.T.
FORMULATRIX®

or



mosquito® LV
SPT Labtech

or



Microlab STAR
Hamilton



epMotion® 5075lc
Eppendorf®

Compared to Olink Explore 3072



Olink®
Accelerating proteomics together

Flexible access to high-quality data via kits and services





Olink®
Accelerating proteomics together

Revealing the truth of human
disease. Protein by protein.



Olink®
Accelerating proteomics together

Protein profiling of ovarian cancer with Olink Explore HT

Background

Unmet need: Accurate early diagnosis of ovarian cancer (OvCa)

Study design: Profile plasma proteins in women with OvCa or benign ovarian tumors using Olink Explore HT and its predecessor, Olink Explore 3072

Study objectives:

- 1) Validate Olink Explore HT's performance by comparing it to Olink Explore 3072
- 2) Investigate the potential of Olink Explore HT to identify a novel OvCa protein signature for more precise diagnostic tests



Ulf Gyllensten, PhD

Professor
Dept of Immunology, Genetics, & Pathology
Uppsala University, Sweden



Olink®
Accelerating proteomics together


Expanded proteome coverage without compromising performance

Technical performance


Compared data from ~2,800 overlapping assays between Olink Explore HT and Olink Explore 3072

Achieved expanded proteome coverage while retaining high data quality

Correlation of Normalized Protein eXpression units (NPX)



Correlation of Δ NPX (effect size) between benign & OvCa groups



Correlation across ~2800 overlapping assays measured by Olink Explore HT and Olink Explore 3072. In the right figure, each data point represents a protein biomarker.




Olink®
Accelerating proteomics together

Expanded coverage yields a significant enrichment of OvCa biomarkers

OvCa-associated proteins

485 proteins biomarkers identified
(Discovery Cohort, n=233)

128 biomarkers replicated (Validation
Cohort, n =171)



OvCa biomarkers identified in the discovery
cohort with Olink Explore platforms




Olink®
Accelerating proteomics together

Increased confidence in biological insights with deeper pathway coverage

Pathway analysis

Higher statistical significance obtained for multiple biological pathways*, providing greater confidence in disease insights.

* Reactome database, <https://reactome.org>



Snippet from gene set enrichment analysis.

Color indicates direction of pathway alterations; color intensity indicates strength of alterations; size of each tile indicates statistical significance



Revealed a novel protein signature that outperformed clinical biomarker

Diagnostic model

Preliminary 19-protein model identified with Olink Explore HT discriminated OvCa from benign samples better than the clinical marker MUC16

3 of these proteins were newly added to the Olink library

<u>Prediction</u>	<u>Signature</u>	<u>Correctly diagnosed</u>
Benign	MUC16	79.1% (68/86)
	19-protein	89.5% (77/86)
OvCa	MUC16	81.2% (69/85)
	19-protein	84.7% (72/85)



Olink®
Accelerating proteomics together

Olink Explore HT: Biomarker discovery with high data quality, throughput, & coverage


Conclusions

High-quality data, yielding similar results as Olink Explore 3072

Expanded proteome coverage, enabling the discovery of additional OvCa biomarkers

Deeper pathway coverage, enhancing statistical confidence in pathway analysis

Superior multi-protein signature than current biomarker, paving the way for improved early OvCa diagnosis





Olink®
Accelerating proteomics together

Protein profiling of cell lysate using Olink Explore HT

Study objectives

1. Evaluate data correlation between Olink Explore HT and its predecessor, Olink Explore 3072
2. Verify the biological relevance of data acquired with Olink Explore HT by comparing protein profiles across different tissues and cancers

Rationale

Cell culture is an excellent model to study biology

Samples

164 cell culture lysate samples: 45 cancers & 22 tissues *



Olink®
Accelerating proteomics together


Excellent correlation across Olink Explore platforms with cell lysate

Technical performance

Compared data between Olink Explore HT and Olink Explore 3072

Analyzed ~2,800 overlapping assays across 28 samples at 3 dilutions (1:1, 1:10, 1:100)

Median correlation coefficient (R) = 0.96



Correlation of overlapping protein assays across all sample dilutions of cell lysate measured with Olink Explore HT and Olink Explore 3072



Olink®
Accelerating proteomics together


Cellular protein profiles accurately reflect tissue- & disease-specific biology (1/2)

UMAP analysis of tissue proteomes

Protein profiles form tissue-specific clusters

Different brain tumors clustered separately

Lymphoid tumors clustered separately from other tissue types






Olink®
Accelerating proteomics together

Cellular protein profiles accurately reflect tissue- & disease-specific biology (2/2)

Hierarchical clustering

6 tissues represented with at least 10 samples included in analysis

Distinct clustering of protein profiles demonstrates tissue and disease discrimination at the cellular level



Hierarchical clustering heatmap of protein profiles from cell cultures representing different tissues and cancers. NPX = Normalized Protein eXpression




Olink®
Accelerating proteomics together

Enriched pathways represent biology of the tissue-of-origin

Pathway analysis

Cellular protein profiles reflect tissue-specific biological pathways

Example: Brain tissue is enriched in neuronal system pathways (yellow), whereas the lymphoid is enriched in immune system pathways (blue)



Enriched pathways for tissue-enriched assays. Color intensity indicates strength of alterations; size of each tile indicates statistical significance.



Olink®
Accelerating proteomics together

Olink Explore HT: A powerful tool to reveal true biological insights

Conclusions

Ensures uncompromised performance, maintaining the same high-quality data as its predecessor

Offers versatile capabilities, enabling high-throughput protein profiling across a wide range of sample types

Facilitates accurate protein measurements, promoting a deeper biological understanding to advance biomarker discovery and drug development